

# SHARE YOUR CGM DATA WITH YOUR CLINIC

DEXCOM CLARITY® CAPTURES GLUCOSE DATA SO YOU AND YOUR CLINIC CAN VIEW YOUR PATTERNS, TRENDS AND STATISTICS ANYTIME, ANYWHERE.

Tech support 1-844-832-1809

## Receive an invitation to share your data.

Your clinic will print or email an invitation to you.

- 1 Click <https://clarity.dexcom.eu/share> or enter it into your internet browser. Confirm to share before the expiration date.
- 2 Click [Log in](#).
- 3 Enter your username and password, then click [Sign in](#).  
Log in with your Dexcom or G5 Mobile account.
- 4 Enter the sharing code from your invitation.
- 5 Enter your date of birth.
- 6 Click [Continue](#).
- 7 Click [I consent to share my data with my clinic](#).
- 8 Click [Yes, Share My Data](#).

**Share data with your clinic**

Sophia Rossi Diabetes Clinic invites you to share your data using Dexcom CLARITY®

Dexcom CLARITY software captures your continuous glucose monitoring (CGM) data so you and your clinic can view patterns, trends and statistics anytime, anywhere.

**Your sharing code\***

**QMCT-PTSX-MMCN**

1 Enter this code at <https://clarity.dexcom.eu/share>

\*Expires: 10/05/2018

**Share Data with a New Clinic**

Your clinic will provide a sharing code for your data so that they have continuous access and can make your appointments more efficient. If your clinic hasn't provided you an invitation to share, you can ask your clinic to create one for you. You can stop sharing with your clinic at any time.

Enter the sharing code provided by your clinic

4

Verify your date of birth

5

6 [Continue](#) [Cancel](#)

**Share Data with a New Clinic**

Please confirm your clinic's name below.

7  I consent to share my data with my clinic.

8 [Yes, Share My Data](#) [Cancel](#)

## To stop sharing with a clinic:

- 1 Log in to CLARITY at [clarity.dexcom.eu](https://clarity.dexcom.eu) with your internet browser.
- 2 Click [Settings](#) at the top of the page.
- 3 Click [Stop Sharing with this Clinic](#) at the bottom of the page.

Data Sharing with Clinics [Share Data with a New Clinic](#)

Access Status	Sharing Data
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

3 [Stop Sharing with this Clinic](#)

**Dexcom**  
**CLARITY®**

LBL015096 Rev002

### BRIEF SAFETY STATEMENT

The web-based Dexcom CLARITY® software is intended for use by healthcare professionals to assist them in the review, analysis, and evaluation of historical CGM data to support effective diabetes management in their patients. It is intended for use as an accessory to Dexcom CGM devices with data interface capabilities. Caution: The software does not provide any medical advice and should not be used for that purpose. Healthcare professionals should use information in the software in conjunction with other clinical information available to them.

© 2017 Dexcom Canada, Inc. All rights reserved.